

# Health-Up T.E.N.S.

## ELECTRONIC ACUPUNCTURE



☯ *Energy from our internal organs is said to circulate around the body through channels called meridians. Imbalances of energy within these meridians indicate illness. Stimulating selected acupuncture points on these meridians can restore the balance of energy by needles or electrical pulses to promote a return to health.*



### FEATURES:

#### ☺ Clinically-proven formula:

By applying to the right acupoints, it generates about 35% effectiveness than that of traditional Acupuncture treatment.

#### ☺ Wireless & Handy & Portable:

Indispensable office-/travel-mate for wellness.

#### ☺ Deep Tissue Acu-"Pressure" through gel:


With water based gel or lotion, the function turns into forceful electrical Acupressure which typical TENS cannot offer.

### MEDICAL-APPLICATION PROPERTIES: FDA 510K: K014126

- ❖ A fixed *1.5Hz* low frequency TENS treatment.
- ❖ Non-invasive, drug-free method of *pain relief*.
- ❖ Symptomatic relief of chronic, intractable pain.
- ❖ Alleviation of muscular pains and neuralgia.



### Unique 2-wave Acu-Pulse



### SPECIFICATIONS:

- Model: MP737-25
- Main device: 1.5 Hz Acu-pulse/TENS Operator, rubber-cap (Rub-1300) inductors.
- Power source: 1 DC 9V (included) battery operated,
- 1 acupoint application chart,
- Carrying leather bag, conduction gel.



### WARNINGS:

1. Qualified physicians' instruction is required to obtain effective treatment.
2. Do not apply the device under the following conditions: (consult qualified physicians as necessary)
  - ☒ Infants or children ☒ Skin disease ☒ Heart disease ☒ Pregnancy ☒ Carrying cardiac pace makers, etc.



ISO 9002  
CE-MDD  
FDA K014126



GMP-Taiwan

### ROLOR ELECTRONICS CORP.

No. 274, Sec. 5, Nanking E. Rd., Taipei, Taiwan  
TEL: 886-2-27682174~5, FAX: 2-27605922~3

E-Mail: webmaster@rolor.com, Web: www.rolor.com